



No further measureable improvements were made with these changes, but these changes may reduce laser jitter etc. Four hoses indicated by colored arrows were moved to burred connectors installed in the capped hole indicated by the same color arrow. The bottom picture shows the final configuration.

